

SAMPLE MASTER Semi-Automatic Biopsy System

Single click for 20mm sample size


Available with or without coaxial introducer needle.

One-handed operation for convenient manipulation of guidance device (ultrasound transducer).

- 20mm Sample Notch
- Echogenic tip for accurate placement under ultrasound guidance
- Centimeter depth markings assist in needle placement
- Can be used with echogenic co-axial introducer needles
- One-handed operation for convenient manipulation of guidance device (ultrasound transducer).

SAMPLE MASTER Disposable Semi-Automatic Biopsy Needle (Box of 10 pcs)

ORDER NO.	DESCRIPTION
SMA 1410	14G x 10cm
SMA 1416	14G x 16cm
SMA 1420	14G x 20cm
SMA 1610	16G x 10cm
SMA 1616	16G x 16cm
SMA 1620	16G x 20cm
SMA 1810	18G x 10cm
SMA 1816	18G x 16cm
SMA 1820	18G x 20cm
SMA 1825	18G x 25cm
SMA 2010	20G x 10cm
SMA 2016	20G x 16cm
SMA 2020	20G x 20cm



20mm Sample Notch

SAMPLE MASTER Co-axial Biopsy System

Allows multiple sampling from one percutaneous puncture with minimal damage to surrounding tissue.

- Trocar point tip maintains directional stability.
- Blunt tip cannula helps prevent inadvertent distal drift and may reduce risk of tissue trauma.
- Echogenic markings facilitate precise ultrasound positioning.
- Locking skin stop helps prevent forward migration.
- Matched in length to cutting cannula.

Co-axial Biopsy Needle for SAMPLE MASTER (Box of 10 pcs)

ORDER NO.	DESCRIPTION
CPX 1410	13G x 5cm, co-axial to SMA 1410
CPX 1416	13G x 11cm, co-axial to SMA 1416
CPX 1420	13G x 15cm, co-axial to SMA 1420
CPX 1610	15G x 5cm, co-axial to SMA 1610
CPX 1616	15G x 11cm, co-axial to SMA 1616
CPX 1620	15G x 15cm, co-axial to SMA 1620
CPX 1810	17G x 5cm, co-axial to SMA 1810
CPX 1816	17G x 11cm, co-axial to SMA 1816
CPX 1820	17G x 15cm, co-axial to SMA 1820
CPX 2010	19G x 5cm, co-axial to SMA 2010
CPX 2016	19G x 11cm, co-axial to SMA 2016
CPX 2020	19G x 15cm, co-axial to SMA 2020



SAMPLE MASTER With Co-axial (Box of 10 pcs)

ORDER NO.	DESCRIPTION
SMAC 1410	SMA 1410 + CPX 1410
SMAC 1416	SMA 1416 + CPX 1416
SMAC 1420	SMA 1420 + CPX 1420
SMAC 1610	SMA 1610 + CPX 1610
SMAC 1616	SMA 1616 + CPX 1616
SMAC 1620	SMA 1620 + CPX 1620
SMAC 1810	SMA 1810 + CPX 1810
SMAC 1816	SMA 1816 + CPX 1816
SMAC 1820	SMA 1820 + CPX 1820
SMAC 2010	SMA 2010 + CPX 2010
SMAC 2016	SMA 2016 + CPX 2016
SMAC 2020	SMA 2020 + CPX 2020



Manufactured by
Adria srl
Via Modena, 46
40017 San Giovanni in Persiceto (BO)
Italy
Tel: +39 0376 559 772

Distributed by
GTA (UK) Ltd
34 Nottingham South Ind. Estate
Ruddington Lane
Wilford
Nottingham, NG11 7EP
Tel: 0115 981 5703
Email: mail@gtamedical.co.uk
Web: www.gta-uk.co.uk