

# HERNIATOME

## Percutaneous Discectomy Device

The Herniectomy, performed with CT and Fluoroscopy or fluoroscopy alone is a minimally invasive procedure on the spine, that aim of which is identical to that of surgery to extract the hernia or a portion of it to reduce the pressure on the nerve root.

- Evaluation in real time of the decompression
- Qualitative and quantitative evaluation
- Cannula with lateral window, that allow a better decompression and extraction of the hernia by covering a wider area with minimal movement
- Available with Curved and Straight Tip
- Available with 20ga introducer needle that ensures the safe and correct positioning of the HERNIATOME

### HERNIATOME Curved Tip (Box of 1)

| ORDER NO. | DESCRIPTION   |
|-----------|---|
| 5091709   | 17G x 9cm, Cervical Type with with 20ga-15cm Introducer |
| 5091715   | 17G x 15cm, Lumbar Type with with 20ga-30cm Introducer  |
| 5091809   | 18G x 9cm, Cervical Type with with 20ga-15cm Introducer |

#### HERNIATOME with curved cannula



#### INTRODUCER



Curved Tip


### HERNIATOME Straight Tip - Direct Access Approach (Box of 1)

| ORDER NO. | DESCRIPTION              |
|-----------|--------------------------|
| 5091749   | 17G x 9cm, Cervical Type |
| 5091745   | 17G x 15cm, Lumbar Type  |
| 5091849   | 18G x 9cm, Cervical Type |

#### HERNIATOME with coclea



#### Cannula with straight tip and stylet



Straight Tip

Manufactured by  
Adria srl  
Via Modena, 46  
40017 San Giovanni in Persiceto (BO)  
Italy  
Tel: +39 0376 559 772

Distributed by  
GTA (UK) Ltd  
34 Nottingham South Ind. Estate  
Ruddington Lane  
Wilford  
Nottingham, NG11 7EP  
Tel: 0115 981 5703  
Email: [mail@gtamedical.co.uk](mailto:mail@gtamedical.co.uk)  
Web: [www.gta-uk.co.uk](http://www.gta-uk.co.uk)