

ISG


Hysterosalpingography Catheter

ISG facilitates injection of contrast media into the uterine cavity

To evaluate patency of the fallopian tubes.

The catheter design allows good physician control while maintaining Sufficient flexibility for comfort.

- The rounded catheter tip and large side holed assures safe placement and superior flow
- Choice of two sizes 5Fr and 7Fr
- Balloon latex free
- The catheter is supplied with introducer, stopcock and syringe



ISG Hysterosalpingography catheter (Box of 5)

| ORDER NO. | DESCRIPTION |
|-----------|-------------|
| 76005 | 5Fr x 30cm |
| 76007 | 7Fr x 30cm |

Manufactured by
Adria srl
Via Modena, 46
40017 San Giovanni in Persiceto (BO)
Italy
Tel: +39 0376 559 772

Distributed by
GTA (UK) Ltd
34 Nottingham South Ind. Estate
Ruddington Lane
Wilford
Nottingham, NG11 7EP
Tel: 0115 981 5703
Email: mail@gtamedical.co.uk
Web: www.gta-uk.co.uk