


RETRIEVE Jamshidi Bone Marrow Biopsy System

RETRIEVE Bone marrow biopsy system are designed to provide quality biopsies following the classical Jamshidi procedure. It is also available in the explant/transplant version.

- Cannula with 5 cutting edges to facilitate penetration
- Stylet with trocar point for smoother entry into the marrow
- Luer lock hub
- Provided with a special cap for a safe removal of the sample
- Available with Explant Version: 5 lateral holes on the external cannula and additional stylet (Code RDSE)



Code RDS



Code RDSE Tip
With 5 side holes


RETRIEVE Bone Marrow Biopsy (Box of 5)

ORDER NO.	DESCRIPTION
RDS 0810	8G x 10cm
RDS 0815	8G x 15cm
RDS 0910	9G x 10cm
RDS 0915	9G x 15cm
RDS 1110	11G x 10cm
RDS 1115	11G x 15cm
RDS 1307	13G x 7cm
RDS 1310	13G x 10cm
RDSE 1110	11G x 10cm, Explant/transplant version
RDSE 1310	13G x 10cm, Explant/transplant version

EVOTRAP Bone Marrow Biopsy System

EVOTRAP Bone marrow biopsy system are designed to provide quality biopsies by Entrapping bone marrow in the needle without altering the sample

- Cannula with 5 cutting edges to facilitate penetration
- Stylet with trocar point for smoother entry into the marrow
- Luer lock hub
- Innovative design provides an exceptional ability to capture and retain bone marrow specimens which eliminates painful deflection technique to harvest a diagnostic quality sample



EVOTRAP System Tip
To retain bone marrow

EVOTRAP Bone Marrow Biopsy (Box of 5)

ORDER NO.	DESCRIPTION
RDT 0910	9G x 10cm
RDT 0915	9G x 15cm
RDT 1110	11G x 10cm
RDT 1115	11G x 15cm
RDT 1310	13G x 10cm

Manufactured by
Adria srl
Via Modena, 46
40017 San Giovanni in Persiceto (BO)
Italy
Tel: +39 0376 559 772

Distributed by
GTA (UK) Ltd
34 Nottingham South Ind. Estate
Ruddington Lane
Wilford
Nottingham, NG11 7EP
Tel: 0115 981 5703
Email: mail@gtamedical.co.uk
Web: www.gta-uk.co.uk